
F · S · T

EN

DE

CN

ES

FR

## REPAIR TIPS FOR TWEEZERS AND FORCEPS

FINE SCIENCE TOOLS

# REPAIR TIPS

## Item No. 29000-00 – Repair Kit for Fine Forceps

### Contains:

- 1 sharpening stone (29008-22) **A**
- 10 O-rings (11200-00) **B**
- 1 bottle instrument oil (50 ml) (29055-00) **C**
- 1 flat nose pliers **D**


Our forceps are highly-precise instruments that have been hand-crafted to meet your needs by experienced surgical technicians – so repair work may require a little practice.

### **We recommend the following procedure:**

Try out the following steps for the first time on old forceps that you don't need anymore.


## STRAIGHTENING FORCEPS

- Using pliers, grip the forceps right at the spot where they are bent.
- Apply more pressure to the tip on the inside (the surface that is touching the other tip) and gently squeeze the pliers.
- Please note that stainless steel and carbon steel are harder alloys and are more likely to break at the bent points than become straight again. Dumoxel, Dumostar and titanium are more flexible alloys.


## ALIGNING FORCEPS

- First, determine where the points of the forceps diverge.
- Then use the same technique on the pliers as in the straightening instructions above, but with more pressure on the side where they are bent – exactly opposite the tip.
- Bend the tip further inwards. The forceps will now close properly again.
- If the tips are not the same length, keep the forceps closed and run the tips gently over a sharpening stone until the protruding end has been removed.


## SHARPENING FORCEPS

- Use the grooved sharpening stone (FST Item No. 29008-22) for sharpening forceps.
- Use an O-ring to keep the tips of the forceps together.
- Hold the forceps gently by the edges and move them gently back and forth on the sharpening stone.
- Always use the instrument oil (FST Item No. 29055-00) when using the sharpening stone.
- Continuously observe the progress of your work under a microscope – and make sure that you sharpen both sides of the edge.


## SHARPENING FORCEPS

- Keep the forceps as parallel to the sharpening stone as possible and repeat the sharpening process on both sides – the more parallel you work, the less abrasion there will be.
- Continuously change the side of the sharpening stone you are using, and clean the dust from the stone with instrument oil every so often.
- Use the grooves in the sharpening stone for the final fine sharpening. To maintain equal tips, rotate the instrument back and forth as you move it over the stone.

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