F · S · T

EN

DE

CN

ES

FR

INSTRUCTION FOR USE OF REUSABLE SURGICAL INSTRUMENTS (IFU)

FINE SCIENCE TOOLS

1. PREPARATION BEFORE USE:

- Thoroughly clean the surgical instruments immediately after each use to prevent the buildup of biofilm and organic material.
- Use an enzymatic detergent solution or an appropriate cleaning agent recommended by the instrument manufacturer.
- Rinse the instruments with distilled or sterile water to remove any remaining detergent residue.
- Inspect the instruments for any visible damage or signs of wear. Do not use instruments that are damaged or compromised in any way.

2. STERILIZATION PROCESS:

- Choose an appropriate sterilization method based on the type of instruments and manufacturer recommendations (e.g. autoclaving, ethylene oxide sterilization).
- Follow the sterilization parameters provided by the instrument manufacturer, including time, temperature, and pressure settings.
- Ensure proper packaging of sterilized instruments to maintain sterility until use. Use sterilization pouches or wraps designed for surgical instruments.

3. STORAGE:

- Store sterilized instruments in a clean, dry, and well-ventilated area.
- Avoid overcrowding or stacking instruments, which may cause damage.
- Consider using instrument trays or containers designed for storage to maintain organization and prevent contamination.

4. PRE-OPERATIVE INSPECTION:

- Before each surgical procedure, inspect the instruments for cleanliness and integrity.
- Check for any residual soil or debris and ensure that the instruments are free from visible contamination.
- Verify that all movable parts are functioning correctly and that there are no signs of corrosion or damage.

5. DURING SURGERY:

- Handle instruments with care to avoid dropping or damaging them.
- Use appropriate techniques for instrument manipulation to minimize the risk of injury to patients or staff.
- Avoid prolonged exposure of instruments to bodily fluids to prevent corrosion or staining.

6. POST-OPERATIVE CARE:

- After use, immediately rinse the instruments with distilled or sterile water to remove any residual blood or tissue.
- Clean the instruments using an enzymatic detergent solution or as per manufacturer instructions to ensure thorough cleaning.
- Inspect the instruments for any damage or signs of wear and report any issues to the appropriate personnel for repair or replacement.

7. REPROCESSING:

- Follow established protocols for reprocessing reusable surgical instruments, including cleaning, sterilization and storage procedures.
- Document each step of the reprocessing process for quality assurance and traceability purposes.
- Regularly review and update reprocessing protocols based on best practices and recommendations from regulatory agencies or professional organizations.

8. MAINTENANCE AND QUALITY CONTROL:

- Implement a regular maintenance schedule for surgical instruments, including sharpening blades and replacing worn components.
- Conduct periodic quality control checks to ensure that instruments meet performance standards and are free from defects.
- Keep detailed records of maintenance activities, including dates of service, repairs, and replacements.

9. TRAINING AND EDUCATION:

- Provide comprehensive training to healthcare personnel involved in the handling, cleaning, and sterilization of surgical instruments.
- Emphasize the importance of following established protocols and guidelines to ensure patient safety and prevent healthcare-associated infections.
- Encourage ongoing education and professional development to stay updated on advances in instrument reprocessing and infection control practices.

10. REPORTING AND FEEDBACK:

- Encourage staff to report any issues or concerns related to the performance or condition of reusable surgical instruments.
- Establish a system for collecting feedback from users to identify areas for improvement in instrument design, reprocessing protocols, or training procedures.
- Investigate and address any reported incidents or adverse events promptly to prevent recurrence and enhance patient safety.

Following these instructions for the proper handling, cleaning, sterilization, and maintenance of reusable surgical instruments will help ensure their safe and effective use in clinical settings, ultimately contributing to improved patient outcomes and satisfaction.

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