


F · S · T

INSTRUMENT *kits*



Our wide range of instruments organized into kits to ensure adequately serving all areas of application

FINE SCIENCE TOOLS

BASIC SURGICAL KIT

ITEM NO. 70006-01

Suitable for common surgeries and tissue handling; a professional-grade set, not a student kit.	1 x	Fine Graefe Forceps – Serrated, Slightly Curved	11151-10
	1 x	Fine Graefe Forceps – 1 x 2 Teeth, Slightly Curved	11154-10
	1 x	Dumont Medical Forceps – Gently Blunted Tips	11253-20
	1 x	Round Handled Micro Needle Holder – Straight, Smooth, 12,5cm, Without Lock	12076-12
	1 x	Ultra Fine Hemostat – Smooth, Curved, 12cm	13021-12
	1 x	Bonn Scissors – Fine Iris Scissors	14084-08
	1 x	Iris Scissors – Adjustable Loops	14106-09
	1 x	Vannas Spring Scissors – 4mm Cutting Edge	15018-10
	1 x	Plastic Sterilization Case made from Ultem – Single Tier	20840-03

ADVANCED SURGICAL KIT

ITEM NO. 70006-02

Intended for users focused on microsurgery and a wider range of advanced procedures, with additional instruments that support easier, more precise tissue handling.	1 x	Grooved Director With Tongue – 13cm	10290-13
	1 x	Fine Graefe Forceps – Smooth, Curved, 10cm	11160-10
	1 x	Dumont Medical Forceps – Gently Blunted Tips	11253-25
	1 x	Dumont #7 Medical Forceps – Tips Close Flat	11274-20
	1 x	Round Handled Micro Needle Holder with Suture Cutter	12075-14
	1 x	Micro-Mosquito Hemostat – Serrated, Straight, 12.5cm	13010-12
	1 x	Spencer Ligature Scissors – For Cutting Stitches	14076-09
	1 x	Bonn Scissors – Highly Polished	14184-09
	1 x	Vannas Spring Scissors – 3mm Cutting Edge	15000-00
	1 x	Goldstein Retractor – 3x3 Sharp, 3cm, 1.8cm Max Spread	17003-03
	1 x	Schwartz Micro Serrefine – Sharp Bend	18052-03
	1 x	Plastic Sterilization Case made from Ultem – Single Tier	20840-03

VASCULAR CATHETERIZATION SURGICAL PACK

ITEM NO. 70003-01

Instrument kit for vascular catheterization in research models up to 100g.	1 x	Grooved Director with tongue	10290-13
	2 x	Graefe Forceps – Serrated, Curved	11051-10
	1 x	Fine Graefe Forceps – 1 × 2 Teeth, Slightly Curved	11154-10
	1 x	Dumont 545 Forceps – Angled 45°	11253-25
	2 x	Moria MC31B Iris Forceps – Smooth, Curved	11373-12
	1 x	Kalt Needle Holder – straight tip	12080-14
	1 x	Micro-Mosquito Hemostat – serrated	13010-12
	1 x	Fine Scissors – Straight	14060-09
	1 x	Fine Scissors – Extra Large Loops	14106-09
	1 x	Vannas-Tübingen Spring Scissors – Sharp angled up tips	15005-08
	1 x	Silk Suture – Size 6/0	18020-60
	1 x	Schwartz Micro Serrefine – light bend	18052-02
	1 x	Absorption Spears – strong, soft, lint-free, biocompatible-proven material	18105-01
	1 x	Magnetic Fixator Retraction System – Elastomer	18200-07
	1 x	Magnetic Fixator Retraction System – Blunt Retractors	18200-10
	1 x	Plastic Sterilization Case – Single Tier	20840-03
	1 x	Assorted Instrument Tip Protectors – vented to prevent rusting	29060-00

VASCULAR CATHETERIZATION SURGICAL PACK

ITEM NO. 70003-02

Instrument kit for vascular catheterization in research models weighing between 100–400g.	2 x	S&T Vessel Dilating Forceps – Angled 45°, 13.5cm, 0.2mm Tip Diameter	00276-13
	1 x	Grooved Director with Tongue – 13cm	10290-13
	1 x	Graefe Forceps – 1x2 Teeth, Curved, 10cm	11054-10
	1 x	Dumont 545 Forceps – Micro-Blunted Standard Tips, Angled 45°, Inox, 11cm	11253-25
	1 x	Kalt Needle Holder – Straight, 13.5cm, Cross-Serrated and Serrated	12080-14
	1 x	Micro-Mosquito Hemostat – Serrated, Straight, 12cm	13010-12
	1 x	Wagner Scissors – Straight, Sharp-Sharp, 12cm	14068-12
	1 x	Strabismus Scissors – Straight, Blunt-Blunt, 11.5cm	14074-11
	1 x	Vannas-Tübingen Spring Scissors – Angled Up, Sharp, 9.5cm, 5mm Cutting Edge	15005-08
	1 x	Suture Thread – Braided Silk, Size 6/0	18020-60
	1 x	Schwartz Micro Serrefine – Light Bend, 26mm	18052-02
	1 x	Absorption Spears – Mounted, Sterile	18105-01
	1 x	Magnetic Fixator Retraction System – Elastomer	18200-07
	1 x	Magnetic Fixator Retraction System – Blunt Retractors, 5mm Wide	18200-11
	1 x	Plastic Sterilization Container – Single Tier	20840-03
	1 x	Instrument Tip Protectors – Assorted Pack	29060-00

RRSSC MICROSURGERY KIT

ITEM NO. 70001-01

A precision instrument kit designed to support delicate microsurgical procedures, based on a cooperation between René Remie Surgical Skill Centre and FST.	1 x	Semken Surgical Forceps – Slim Shape, 1x2 Teeth	11024-12
	1 x	Graefe Forceps – 1x2 Teeth, Straight, 10cm	11053-10
	1 x	Extra Fine Graefe Forceps – Serrated, Straight, 10cm	11150-10
	1 x	Extra Fine Graefe Forceps – Serrated, Curved, 10cm	11152-10
	1 x	Dumont 5 Forceps – Micro-Blunted Standard Tips, Straight, Inox, 11cm	11253-20
	1 x	Dumont 545 Forceps – Micro-Blunted Standard Tips, Angled 45°, Inox, 11cm	11253-25
	1 x	Round Handled Needle Holder – Curved Tips	12074-12
	1 x	Crile-Wood Needle Holder – Tungsten Carbide, Straight, Serrated, 15cm, with Lock	12503-15
	1 x	Hartmann Hemostat – Serrated, Straight, 10cm	13002-10
	1 x	Surgical Scissors – Tungsten Carbide, Straight, Sharp-Sharp, 14.5cm	14502-14
	1 x	Fine Scissors – Tungsten Carbide, Straight, Sharp-Sharp, 11.5cm	14568-12
	1 x	Vannas-Tübingen Spring Scissors – Angled Up, Sharp, 9.5cm, 5mm Cutting Edge	15005-08
	1 x	Spring Scissors – ToughCut®, Straight, Sharp, 12.5cm, 6mm Cutting Edge	15124-12
	1 x	Bulldog Serrefine – Straight, 35mm	18050-35
	1 x	Biemer Clip – Straight, 6mm	18058-15
	1 x	Biemer Clip – Straight, 9mm	18058-18
	1 x	Biemer Clip Applicator Forceps Style – Straight, without Lock, 14 cm	18059-14
	1 x	Vessel Cannulation Forceps – 0.35mm OD, 11.5cm	18403-11
	1 x	Sterilization Container – Double Tier (with RRSSC logo)	24842-03

PROFESSIONAL ZEBRAFISH KIT

ITEM NO. 70004-01

For advanced users, with essential tools for common Zebrafish procedures.	1 x	Probe – Straight, 15.5cm, 0.25mm Tip Diameter	10140-01
	1 x	Probe – Angled 80° Short, 15.5cm, 0.15mm Tip Diameter	10140-02
	1 x	Tissue Forceps – 1x2 Teeth, Straight, 12cm	11021-12
	1 x	Graefe Forceps – Serrated, Curved, 10cm	11051-10
	1 x	Dumont #5 Forceps – Standard Tips, Straight, Inox, 11cm	11251-20
	1 x	Noyes Spring Scissors – Straight, Sharp, 12cm, 14mm Cutting Edge	15012-12
	1 x	Plastic Sterilization Container - Single Tier	20840-02

BASIC ENTOMOLOGY KIT ITEM NO. 70007-01


Instrument kit for systematic and morphological studies, featuring narrow insect-handling forceps for most common specimens.

1 x	Probe – Straight, 15.5cm, 0.25mm Tip Diameter	10140-01
1 x	Kuehne Forceps – Straight, Smooth	11074-02
1 x	Dumont #5SF Forceps – Straight, Super Fine	11252-00
1 x	Dumont #7 Medical Forceps – Curved	11274-20
1 x	Colored Plastic Instrument Case – Single Tier	20830-00
1 x	Insect Pins – Black Anodized	26000-45
1 x	Live Insect Forceps – Narrow, Straight	26029-10

ADVANCED ENTOMOLOGY KIT ITEM NO. 70007-02

Instrument kit for complex research and specimen preparation, with narrow and wide forceps for a broader range of insect sizes and additional tools for dissection.

1 x	Side Cut Micro Knife – Cutting Edge on the Angled Side	10058-15
1 x	Wecker Iris Spatulae – Curved, 12cm	10092-12
1 x	Probe – Straight, 15.5cm, 0.25mm Tip Diameter	10140-01
1 x	Vial Forceps – Straight, Smooth	11006-15
1 x	Graefe Forceps – Slightly Curved, Serrated	11051-10
1 x	Kuehne Forceps – Straight, Smooth	11074-02
1 x	Dumont #5SF Forceps – Straight, Super Fine	11252-00
1 x	Dumont #7 Medical Forceps – Curved	11274-20
1 x	#5 Self-Closing Fine Forceps – Standard Tips, Straight, 11cm	11485-11
1 x	Vannas Spring Scissors – 2mm Cutting Edge	15000-03
1 x	Plastic Sterilization Case made from Ultem – Single Tier	20840-03
1 x	Insect Pins – Black Anodized	26000-45
1 x	Insect Pins – Black Anodized	26000-50
1 x	Live Insect Forceps – Narrow, Straight	26029-10
1 x	Live Insect Forceps – Wide, Straight	26030-10



LOOKING FOR A DIFFERENT COMBINATION?
Custom kits can be tailored to your application.
Contact your FST representative or
authorized reseller to explore the options.

**DESIGN
YOUR OWN
KIT**