



S&T Vessel Dilating Forceps - Angled 10°/11cm/0.2mm Tip Diameter



Designed with elongated, highly polished, parallel tips for controlled intraluminal vessel dilation.

General Info

Product Number	00125-11
Category	Forceps & Tweezers > Smooth
EAN	4251909500061



Find 00125-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Tipdiameter	0.2 mm		
Tipdimensions	0.2 x 0.1 mm		
Tips	High polished, smooth		
Tipshape	Angled 10°		