



S&T Vessel Dilating Forceps- Angled 45°/13.5cm/0.3mm Tip Diameter



Designed with elongated, highly polished, parallel tips for controlled intraluminal vessel dilation.

General Info

Product Number	00276-13
Category	Forceps & Tweezers > Smooth
EAN	4251909500092



Find 00276-13 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13.5 cm		
Tipdiameter	0.2 mm		
Tipdimensions	0.3 x 0.2 mm		
Tips	High polished, smooth		
Tipshape	Angled 45°		