



S&T Vascular Clamp - Straight/7mm



Allgemein

Produktnummer	00396-01
Kategorie	Clamps & Clips > Clamps
EAN	4251909500115



Finde 00396-01 auf finescience.com

Technische Informationen

Versandgewicht	5 g	Nettogewicht	1 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	3.8 mm		
Clamppress	Variable		
Length	7 mm		
Grampress	5 to 15 g		
Jawdims	3.5 x 1 mm		
Serrated	No		
Shape	Straight		