



S&T Vascular Clamp - Straight/7mm



General Info

Product Number	00396-01
Category	Clamps & Clips > Clamps
EAN	4251909500115



Find 00396-01 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	3.8 mm		
Clamppress	Variable		
Length	7 mm		
Grampress	5 to 15 g		
Jawdims	3.5 x 1 mm		
Serrated	No		
Shape	Straight		