



S&T Vascular Clamp - Straight/11mm



Allgemein

Produktnummer	00398-02
Kategorie	Clamps & Clips > Clamps MedTech Line
EAN	4251909500122



Finde 00398-02 auf
finescience.com

Technische Informationen

Versandgewicht	5 g	Nettogewicht	1 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	5.5 mm		
Clamppress	Variable		
Length	11 mm		
Grampress	5 to 15 g		
Jawdims	5.5 x 1.5 mm		
Serrated	No		
Shape	Straight		