



S&T Vascular Clamp - Straight/16mm



General Info

Product Number	00400-03
Category	Clamps & Clips > Clamps
EAN	4251909500139



Find 00400-03 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	7.5 mm		
Clamppress	Variable		
Length	16 mm		
Grampress	5 to 15 g		
Jawdims	7.5 x 1.75 mm		
Serrated	No		
Shape	Straight		