



S&T Vessel Cannulation Forceps - 0.6mm OD/11cm



These vessel cannulation forceps are useful for introducing a small plastic tube into a blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is guided into the vessel opening while preventing damage from unwanted movement.

General Info

Product Number	00574-11
Category	Forceps & Tweezers > Cannulation
EAN	4251909500146



Find 00574-11 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Tipdimensions	1 x 1.2 mm		
Tips	Hollow		
Tipshape	Angled		
Tubingsize	0.6 mm OD		