



S&T Vessel Cannulation Forceps - 0.35 OD/11cm



These vessel cannulation forceps are useful for introducing a small plastic tube into a blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is guided into the vessel opening while preventing damage from unwanted movement.

General Info

|                |                                  |
|----------------|----------------------------------|
| Product Number | 00608-11                         |
| Category       | Forceps & Tweezers > Cannulation |
| EAN            | 4251909508876                    |



Find 00608-11 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |     |
|------------------------|---------------------------------|-----------------|-----|
| Shipping Weight        | 30 g                            | Net Weight      | 9 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA  |
| Minimum Order Quantity | 1                               |                 |     |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |     |
| Length                 | 11 cm                           |                 |     |
| Tipdimensions          | 0.6 x 0.5mm                     |                 |     |
| Tips                   | Hollow                          |                 |     |
| Tipshape               | Angled                          |                 |     |
| Tubingsize             | 0.35 mm OD                      |                 |     |