



SuperGrip tips are processed with MicroCrystals, a method developed by S&T to infuse Tungsten Carbide on stainless steel. The formation of this ultra-fine coating guarantees atraumatic handling of tissue while improving the grip dramatically. The gripping surface is similar to a match head. Recommended when forceps are utilized as needle holders or when slippery tissues might be encountered.

General Info

Product Number	00632-11
Category	Forceps & Tweezers > Smooth
EAN	4251909500153



Find 00632-11 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide tips		
Length	11 cm		
Style	JF-5TC		
Tipdimensions	0.3 x 0.25 mm		
Tips	SuperGrip		
Tipshape	Straight		