

F · S · T

Disposable Safety Scalpels #11 / Angled / 20mm cutting edge



The protective cartridge completely covers the blade for safe, fast, and easy disposal. Eliminates accidental transmission of blood borne pathogens. Designed with a convenient ruler for quick measurement. Single-use and non-autoclavable.

General Info

Product Number	10000-11
Category	Scalpels & Knives > Disposable Scalpels
EAN	787551012200



Find 10000-11 on
finescience.com

Technical Info

Shipping Weight	120 g	Net Weight	60 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with plastic handle		
Autoclavesafe	No		
Bladeshape	Angled		
Cuttingedge	20 mm		
Length	14 cm		
Sterile	No		
Thickness	0.4 mm		