



Disposable Safety Scalpels #15 / Curved / 15mm cutting edge



The single-use protective cartridge completely covers the blade for safe, fast, and easy disposal. Eliminates accidental transmission of blood borne pathogens. Designed with a convenient ruler for quick measurement.

General Info

Product Number	10000-15
Category	Scalpels & Knives > Disposable Scalpels
EAN	787551012217



Find 10000-15 on
finescience.com

Technical Info

Shipping Weight	120 g	Net Weight	60 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with plastic handle		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	15 mm		
Length	14 cm		
Sterile	No		
Thickness	0.4 mm		