



Scalpel Blade Removal & Disposal Box

Used to safely remove and dispose contaminated scalpel blades. Can be used for scalpels that fit the #3, #4, #7 and gross anatomy handles.



General Info

Product Number	10002-50
Category	Scalpels & Knives > Blade Remover
EAN	4251909500184



Find 10002-50 on
finescience.com

Technical Info

Shipping Weight	265 g	Net Weight	223 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Capacity	80 Blade		
Length	7.5 cm		
Width	2.5 cm		
Height	5 cm		
Dimensionstext	7.5 x 5 x 2.5 cm		