



A safety lock holds the blade securely, and allows simple and safe blade attachment and removal.

### General Info

Product Number	10003-13
Category	Scalpels & Knives > Scalpel Handles
EAN	4251909500191



Find 10003-13 on  
[finescience.com](https://www.finescience.com)

### Technical Info

Shipping Weight	30 g	Net Weight	30 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	12 cm		
Use with	Scalpel blades		