



The Safety Scalpel Handle #3 is compatible with all Safety Scalpel Blades. The handle is reusable and made of Teflon coated stainless steel which provides a non-stick surface. The Blade cartridges are color-coded with clear windows for easy identification. The protective cartridge completely covers the blade for safe, fast, and easy disposal. Handle is reusable and autoclavable. Blades are single-use and non-sterile. Handle and cartridges sold separately.

General Info

Product Number	10006-12
Category	Scalpels & Knives > Scalpel Handles
EAN	4251909500214



Find 10006-12 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Teflon coated Surgical Stainless Steel (Inox)		
Length	12.5 cm		