



Scalpel Blades - #12



Blades come individually wrapped in transparent packages. The alloy has a high carbon content which provides a longer lasting cutting edge. Not autoclavable.



General Info

Product Number	10012-00
Category	Scalpels & Knives > Blades
EAN	4061706158446



Find 10012-00 on
finescience.com

Technical Info

Shipping Weight	105 g	Net Weight	92 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	15 mm		
Sterile	No		
Style	#12 blades		
Thickness	0.4 mm		