



Safety Scalpel Blade Cartridges - #15



The Safety Scalpel Handle #3 is compatible with all Safety Scalpel Blades. The handle is reusable and made of Teflon coated stainless steel which provides a non-stick surface. The Blade cartridges are color-coded with clear windows for easy identification. The protective cartridge completely covers the blade for safe, fast, and easy disposal. Handle is reusable and autoclavable. Blades are single-use and non-sterile. Handle and cartridges sold separately.

General Info

Product Number	10015-10
Category	Scalpels & Knives > Blades
EAN	787551012170



Find 10015-10 on finescience.com

Technical Info

Shipping Weight	120 g	Net Weight	20 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	15 mm		
Sterile	No		
Thickness	0.4 mm		