



Blades come individually wrapped in transparent packages. The alloy has a high carbon content which provides a longer lasting cutting edge. Not autoclavable.



General Info

Product Number	10020-00
Category	Scalpels & Knives > Blades
EAN	4251909500252



Find 10020-00 on
finescience.com

Technical Info

Shipping Weight	140 g	Net Weight	121 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	30 mm		
Sterile	No		
Style	#20 blades		
Thickness	0.4 mm		