



The alloy used has a high carbon content which provides a longer lasting edge. Blades come individually wrapped.

General Info

Product Number	10022-00
Category	Scalpels & Knives > Blades
EAN	4251909500276



Find 10022-00 on finescience.com

Technical Info

Shipping Weight	185 g	Net Weight	166 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	40 mm		
Sterile	No		
Style	#22 blades		
Thickness	0.4 mm		