




Ceramic Scalpel Blades - Angled 45°/2.5 mm



These ceramic blades have cutting edges that stay sharp 100 times longer than steel, are rust and corrosion free, acid and chemical resistant, autoclavable, non-magnetic, and without trace metal contamination.

General Info

Product Number	10025-43
Category	Scalpels & Knives > Blades
EAN	4251909512132



Find 10025-43 on
finescience.com

Technical Info

Shipping Weight	4 g	Net Weight	2 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Ceramic		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Bladeshape	Angled 45°		
Cuttingedge	2.5 mm		
Dimensionstext	33 x 6 x 0.5 mm		
Sterile	No		
Thickness	0.5 mm		