




Ceramic Blades - Angled / 8mm Cutting Edge (Pack of 5)



These ceramic blades have cutting edges that stay sharp 100 times longer than steel, are rust and corrosion free, acid and chemical resistant, autoclavable, non-magnetic, and without trace metal contamination.

General Info

Product Number	10025-45
Category	Scalpels & Knives > Blades
EAN	4251909500290



Find 10025-45 on
finescience.com

Technical Info

Shipping Weight	5 g	Net Weight	1 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Ceramic		
Antimagnetic	Yes		
Autoclavesafe	Yes		
Bladeshape	Angled 45°		
Cuttingedge	8 mm		
Dimensionstext	33 x 6 x 0.5 mm		
Sterile	No		
Thickness	0.5 cm		