



Ceramic Blade Cutter



Ceramic edge stays sharp 100 times longer than a steel edge. Rust and corrosion free. No trace metal contamination. The handle is not autoclavable. The blade is autoclavable.

General Info

Product Number	10025-70
Category	Scalpels & Knives > Knives
EAN	4251909500313



Find 10025-70 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	23 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic handle with ceramic blade		
Antimagnetic	Yes		
Autoclavesafe	Partial		
Bladeshape	Straight		
Cuttingedge	58 mm		
Length	13 cm		
Sterile	No		
Thickness	0.3 mm		