



Tübingen Micro Probe - Angled 45°/13.5cm/0.03mm Tip Diameter



General Info

Product Number	10032-13
Category	Pins, Probes & Hooks > Probes
EAN	4251909500344



Find 10032-13 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	40 g	Net Weight	2 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13.5 cm		
Tipdiameter	0.03 mm		
Tipshape	Angled 45°		