



Interchangeable Blades - Angled/10mm Cutting Edge



General Info

Product Number	10035-15
Category	Scalpels & Knives > Blades
EAN	4251909500405



Find 10035-15 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Angled		
Cuttingedge	10 mm		
Sterile	No		
Thickness	0.6 mm		