



Interchangeable Blades - Curved/10mm Cutting Edge



General Info

Product Number	10035-20
Category	Scalpels & Knives > Blades
EAN	4251909500412



Find 10035-20 on
finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	7 g
Package Size	5 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	10 mm		
Sterile	No		
Thickness	0.6 mm		