



Scalpel Blades - Breakable/Variable Cutting Edge



These are the only blades that break cleanly. The alloy has a high carbon content which provides a longer lasting cutting edge. Not autoclavable.



General Info

Product Number	10050-00
Category	Scalpels & Knives > Blades
EAN	4251909500443



Find 10050-00 on [finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Straight		
Bladethickness	0.1 mm		
Cuttingedge	35 mm		
Sterile	No		