



Blade Holder & Breaker - Flat Jaws/11cm



Used for producing custom sized microblades.

General Info

Product Number	10052-11
Category	Scalpels & Knives > Scalpel Handles
EAN	4251909500450



Find 10052-11 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	26 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Jawwidth	4 mm		
Tipshape	Flat		