



Blade Holder & Breaker - Concave-Convex Jaws/9cm



Used for producing custom sized microblades The concave and convex jaws provide a secure hold. Round handles allow easy fingertip adjustments.

General Info

Product Number	10053-09
Category	Scalpels & Knives > Scalpel Handles
EAN	4251909500467



Find 10053-09 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	9 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	9 cm		
Handlediameter	5.5 mm		
Jawwidth	4 mm		
Tipshape	Concave / Convex		