



Blade Holder & Breaker - Concave-Convex Jaws/13cm



Used for producing custom sized microblades The concave and convex jaws provide a secure hold. Round handles allow easy fingertip adjustments.

General Info

Product Number	10053-13
Category	Scalpels & Knives > Scalpel Handles
EAN	4251909500474



Find 10053-13 on
finescience.com

Technical Info

Shipping Weight	60 g	Net Weight	26 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Handlediameter	8 mm		
Jawwidth	4 mm		
Tipshape	Concave / Convex		