



Dissecting Knife - Angled 30°/12.5cm



The cutting edge is lance shaped to provide the ideal cutting angle.

General Info

Product Number	10056-12
Category	Scalpels & Knives > Micro Knives
EAN	4251909500498



Find 10056-12 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladeshape	Angled 30°		
Cuttingedge	4 mm		
Length	12.5 cm		
Thickness	0.5 mm		