



Micro Knife - Side Cut/15cm



The cutting edge is on the angled side.

General Info

Product Number	10058-15
Category	Scalpels & Knives > Micro Knives
EAN	4251909500504



Find 10058-15 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladeshape	Side cut		
Cuttingedge	3 mm		
Length	15 cm		
Thickness	0.2 mm		