



Graefe Knife - Angled 90°/13cm



Ideal for cutting dura. The cutting edge is on the lower side of the tip.

General Info

Product Number	10059-13
Category	Scalpels & Knives > Micro Knives
EAN	4251909500511



Find 10059-13 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladeshape	Angled 90°		
Cuttingedge	1.5 mm		
Length	13 cm		
Thickness	0.5 mm		