



Circle Knife - Round/15cm



The cutting edge wraps completely around the 2mm tip diameter.

General Info

Product Number	10059-15
Category	Scalpels & Knives > Micro Knives
EAN	4251909500528



Find 10059-15 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladeshape	Round		
Cuttingedge	2 mm		
Length	16 cm		
Thickness	0.5 mm		