



These gross anatomy blades come individually wrapped in transparent packages. The alloy has a high carbon content which provides a longer lasting cutting edge. Not autoclavable.

General Info

Product Number	10060-00
Category	Scalpels & Knives > Blades
EAN	4061706144722



Find 10060-00 on
finescience.com

Technical Info

Shipping Weight	310 g	Net Weight	291 g
Package Size	100 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Carbon steel		
Autoclavesafe	No		
Bladeshape	Curved		
Cuttingedge	60 mm		
Length	7.5 cm		
Sterile	No		
Style	#60 blades		
Thickness	0.6 mm		