



Scalpel Handle - #4 Gross Anatomy



This scalpel handle will only hold the #60 gross anatomy blades.

General Info

|                |                                     |
|----------------|-------------------------------------|
| Product Number | 10060-13                            |
| Category       | Scalpels & Knives > Scalpel Handles |
| EAN            | 4061706014261                       |



Find 10060-13 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 30 g                            | Net Weight      | 27 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 13 cm                           |                 |      |