



Micro Probe - Angled 90° Short/15.5cm/0.03mm Tip Diameter



Allgemein

|               |                               |
|---------------|-------------------------------|
| Produktnummer | 10063-15                      |
| Kategorie     | Pins, Probes & Hooks > Probes |
| EAN           | 4251909500559                 |

Finde 10063-15 auf  
[finescience.com](https://finescience.com)

Technische Informationen

|                     |                                 |              |     |
|---------------------|---------------------------------|--------------|-----|
| Versandgewicht      | 40 g                            | Nettogewicht | 4 g |
| Paketgröße          | 1 Pc.                           | Maßeinheit   | EA  |
| Mindestbestellmenge | 1                               |              |     |
| Alloymaterial       | Surgical Stainless Steel (Inox) |              |     |
| Length              | 15.5 cm                         |              |     |
| Tipdiameter         | 0.03 mm                         |              |     |
| Tipshape            | Angled 90° Short                |              |     |