



Sichel Knife - Curved/12cm



General Info

Product Number	10073-14
Category	Scalpels & Knives > Micro Knives
EAN	4251909500597



Find 10073-14 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladeshape	Curved		
Cuttingedge	9 mm		
Length	12 cm		
Thickness	0.6 mm		