



Delicate Bone Scraper - Angled/16.5cm/0.4mm Cutting Edge



The edges of this instrument are medium sharp so that they can scrape tissue from delicate bone without cutting into the bone itself.

General Info

Product Number	10075-16
Category	Chisels & Currettes
EAN	4251909500603



Find 10075-16 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	3.5 mm		
Length	16.5 cm		
Tipshape	Angled		