



Micro Curette - 13cm/1mm Tip Diameter



A curette tip consists of a scoop with sharp edges and is used for scraping or cleaning. It also can be used for separating and transporting, inoculating plates or transferring small amounts of chemicals.

General Info

Product Number	10081-10
Category	Chisels & Curettes
EAN	4251909500627



Find 10081-10 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	13 cm		
Serrated	No		
Tipdiameter	1 mm		