



Micro Curette - 13cm/1.5mm Tip Diameter



A curette tip consists of a scoop with sharp edges and is used for scraping or cleaning. It also can be used for separating and transporting, inoculating plates or transferring small amounts of chemicals.

General Info

Product Number	10082-15
Category	Chisels & Curettes
EAN	4251909500634



Find 10082-15 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	15 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Serrated	No		
Tipdiameter	1.5 mm		