



Micro Curette - 12.5cm/2mm Tip Diameter/Serrated



A curette tip consists of a scoop with sharp edges and is used for scraping or cleaning. It also can be used for separating and transporting, inoculating plates or transferring small amounts of chemicals.

General Info

Product Number	10083-20
Category	Chisels & Curettes
EAN	4251909500641



Find 10083-20 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	2		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Serrated	Yes		
Tipdiameter	2 mm		