



Double Sided Micro Spatula - 13cm



Allgemein

Produktnummer	10091-12
Kategorie	Spoons & Spatulae
EAN	4251909500702



Finde 10091-12 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	13 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	25 x 2 x 0.3 mm		
Tipthickness	0.3 mm		
Tipwidth	2 mm		
Tipshape	Curved/Angled		