



Double Sided Micro Spatula - 13cm



General Info

Product Number	10091-12
Category	Spoons & Spatulae
EAN	4251909500702



Find 10091-12 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	25 x 2 x 0.3 mm		
Tipthickness	0.3 mm		
Tipwidth	2 mm		
Tipshape	Curved/Angled		