



Wecker Iris Spatulae - Slight Curve/13cm



Allgemein

Produktnummer	10093-13
Kategorie	Spoons & Spatulae Most Popular > Spoons & Spatulae
EAN	4251909500726



Finde 10093-13 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	15 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12.5 cm		
Tipdimensions	25 x 1.8 x 0.5 mm		
Tipthickness	0.5 mm		
Tipwidth	1.8 mm		
Tipshape	Slightly curved		