



Wecker Iris Spatulae - Straight/13cm



Allgemein

Produktnummer	10094-13
Kategorie	Spoons & Spatulae Most Popular > Spoons & Spatulae
EAN	4251909500733



Finde 10094-13 auf
finescience.com

Technische Informationen

Versandgewicht	40 g	Nettogewicht	17 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	25 x 3 x 0.5 mm		
Tipthickness	0.5 mm		
Tipwidth	3 mm		
Tipshape	Straight		