



Wecker Iris Spatulae - Straight/13cm



General Info

Product Number	10094-13
Category	Spoons & Spatulae Most Popular > Spoons & Spatulae
EAN	4251909500733



Find 10094-13 on
finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	13 cm		
Tipdimensions	25 x 3 x 0.5 mm		
Tipthickness	0.5 mm		
Tipwidth	3 mm		
Tipshape	Straight		