



Dissecting Bone Chisel - Angled/12.5cm/5mm Cutting Edge



Allgemein

Produktnummer	10095-12
Kategorie	Chisels & Curettes
EAN	4251909500740



Finde 10095-12 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	15 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	12.5 cm		
Tipthickness	0.8 mm		
Tipshape	Angled		