



Austin Bone Chisel - Straight/16cm/1.5mm Cutting Edge



This product is discontinued and has limited availability.

General Info

Product Number	10097-16
Category	Chisels & Curettes
EAN	4251909500764



Find 10097-16 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	1.5 mm		
Length	16 cm		
Tipthickness	0.05 mm		
Tipshape	Straight		